
Rabbit Anti-TNF Antibody

Cat# AB3001-100, AB3001-500

Product description

Product: Rabbit Anti-TNF Antibody
Host: Rabbit
Species: Human
Isotype: IgG
Format: Purified
Immunogen: Recombinant Human TNF Alpha
Product Type: Polyclonal Antibody

Formulation

Rabbit Anti-TNF Antibody stored in 1X Phosphate Buffer Saline, 20% Glycerol.

Application

The product is purified on protein-G affinity resin and specifically recognizes the Human TNF alpha. The antibody can be used in bioanalytical assay and western blot. It is recommended that the end user optimizes the product for their particular application with appropriate negative and positive controls.

Specificity

This Rabbit Anti TNF antibody recognizes specifically the Human TNF Alpha.

Tested Application

This product has been tested for ELISA and western blot application.

Storage

Store at: -20°C in aliquots and avoid repeated freeze/thaw cycles.

Pack size available

100 µg, 500 µg

Regulatory:

For research use only (RUO)

This Product is for Research Use only (RUO) and not to be used in Diagnostic Procedure.