

## Rabbit Anti-Adalimumab Antibody

Cat# AB3004\_100, AB3004\_500

### Product description

---

Product: Rabbit Anti-Adalimumab Antibody  
Format: Purified  
Immunogen: Adalimumab  
Product Type: Polyclonal Antibody  
Other Names: Humira®  
Quantity: 100 & 500 µg

### Formulation

---

Rabbit Anti-Adalimumab Antibody stored in 1X Phosphate Buffer Saline, 20% Glycerol.

### Application

---

The product is purified on protein A affinity resin and recognize Adalimumab in indirect ELISA application. It is recommended that the end user optimizes the product for their particular application with appropriate negative and positive controls.

### Storage

---

Store at: -20°C in aliquots and avoid repeated freeze/thaw cycles.

### Pack size available

---

100 µg, 500 µg

### Regulatory:

---

For research use only (RUO)

*This Product is for Research Use only (RUO) and not to be used in Diagnostic Procedure.*