

## Mouse Anti-TNF Antibody

Cat# AB3005-100, AB3005-500

### Product description

---

Product: Mouse Anti-TNF Antibody  
Host: Mouse  
Species: Human  
Isotype: IgG  
Format: Purified  
Immunogen: Recombinant Human TNF Alpha  
Product Type: Monoclonal Antibody

### Formulation

---

Mouse Anti-TNF Antibody stored in 1X Phosphate Buffer Saline, 20% Glycerol.

### Application

---

The product is purified on protein-G affinity resin and specifically recognizes the Human TNF alpha. The antibody can be used in bioanalytical assay and western blot. It is recommended that the end user optimizes the product for their particular application with appropriate negative and positive controls.

### Specificity

---

This Mouse Anti TNF antibody recognizes specifically the Human TNF Alpha.

### Tested Application

---

This product has been tested for ELISA and western blot application.

### Storage

---

Store at: -20°C in aliquots and avoid repeated freeze/thaw cycles.

### Pack size available

---

100 µg, 500 µg

### Regulatory:

---

For research use only (RUO)

*This Product is for Research Use only (RUO) and not to be used in Diagnostic Procedure.*