

Rabbit Anti-Bevacizumab Antibody

Cat# AB3006-100, AB3006-500

Product description

Product: Rabbit Anti-Bevacizumab Antibody
Format: Purified by Depletion against Human IgG
Immunogen: Bevacizumab (Fab)₂
Product Type: Polyclonal Antibody
Other Names: Anti-Avastin® Antibody

Formulation

Rabbit Anti-Bevacizumab Antibody stored in 1X Phosphate Buffer Saline, 20% Glycerol.

Application

The product is purified on protein-A affinity resin and depleted against Human IgG. This antibody recognizes Bevacizumab in ELISA application.

It is recommended that the end user optimizes the product for their particular application with appropriate negative and positive controls.

Storage

Store at -20°C in aliquots and avoid repeated freeze/thaw cycles.

Pack size available

100 µg, 500 µg

Regulatory:

For research use only (RUO)

This Product is for Research Use only (RUO) and not to be used in Diagnostic Procedure.