

## Mouse Anti-GCSF Antibody

Cat# AB3007-100, AB3007-500

### Product description

---

Product: Mouse Anti-GCSF Antibody  
Format: Purified  
Immunogen: Pegylated GCSF protein  
Product Type: Monoclonal Antibody  
Other Names: Peg-grafeel®

### Formulation

---

Mouse Anti-GCSF Antibody stored in 50mM Phosphate Buffer Saline, 20% Glycerol.

### Application

---

The product is purified on protein G affinity resin and recognize human GCSF in ELISA application. It is recommended that the end user optimizes the product for their particular application with appropriate negative and positive controls.  
ELISA: Use at 1µg/ml Western: Use at 5-10 µg/ml

### Storage

---

Store at: -20°C in aliquots and avoid repeated freeze/thaw cycles.

### Pack size available

---

100 µg, 500 µg

### Regulatory:

---

For research use only (RUO)

*This Product is for Research Use Only (RUO) and not to be used in diagnostic procedure*