

Mouse Anti-PEG Antibody

Cat# AB3008-100, AB3008-500

Product description

Product: Mouse Anti-PEG Antibody
Format: Purified
Immunogen: Pegylated PEG protein
Product Type: Monoclonal Antibody
Other Names: Peg-grafeel®

Formulation

Mouse Anti-PEG Antibody stored in 50mM Phosphate Buffer Saline, 20% Glycerol.

Application

The product is purified on protein G affinity resin and recognize PEG in ELISA application. It is recommended that the end user optimizes the product for their particular application with appropriate negative and positive controls.

Tested Application

ELISA: Use at 1 µg/ml
Western: Use at 5-10 µg/ml

Storage

Store at: -20°C in aliquots and avoid repeated freeze/thaw cycles.

Pack size available

100 µg, 500 µg

Regulatory:

For research use only (RUO)

This Product is for Research Use Only (RUO) and not to be used in diagnostic procedure