

Anti-Ranibizumab Polyclonal Antibody

Cat# AB3012-100, AB3012-500

Product description

Product: Anti-Ranibizumab Polyclonal Antibody
Host: Rabbit
Isotype: IgG
Format: Protein-A affinity Purified and depleted against Human IgG.
Immunogen: Ranibizumab
Product Type: Polyclonal Antibody
Other Names: Lucentis®

Formulation


Anti-Ranibizumab Polyclonal Antibody stored in 1X Phosphate Buffer Saline, 20% Glycerol.

Application

The product is purified on protein-A affinity resin and specifically recognizes the Ranibizumab drug. The antibody can be used as positive control for immunogenicity assay and for western blot. It is recommended that the end user optimizes the product for their particular application with appropriate negative and positive controls.

Specificity

This Anti-Ranibizumab Polyclonal Antibody recognizes Ranibizumab specifically.



Tested Application

This product has been tested for ELISA format and western blot application.

Storage

Store at: -20°C in aliquots and avoid repeated freeze/thaw cycles.

Pack size available

100 µg, 500 µg

Regulatory:

For research use only (RUO)

This Product is for Research Use only (RUO) and not to be used in Diagnostic Procedure.