

Rabbit Anti-Trastuzumab Antibody

Cat# AB3014-100, AB3014-500

Product description

Product: Rabbit Anti-Trastuzumab Antibody

Host: Rabbit

Isotype: IgG

Format: Protein-A affinity Purified and depleted against Human IgG.

Immunogen: Trastuzumab

Product Type: Polyclonal Antibody

Other Names: Herceptin®

Formulation

Rabbit Anti-Trastuzumab Antibody stored in 50mM Phosphate Buffer Saline, 20% Glycerol.

Application

The product is purified on protein-A affinity resin and specifically recognizes the Trastuzumab drug. The antibody can be used as positive control for immunogenicity assay and for western blot. It is recommended that the end user optimizes the product for their particular application with appropriate negative and positive controls.

Specificity

This Rabbit Anti-Trastuzumab Antibody recognizes Trastuzumab specifically.

Tested Application

This product has been tested for ELISA format and western blot application.

Storage

Store at: -20°C in aliquots and avoid repeated freeze/thaw cycles.

Pack size available

100 µg, 500 µg

Regulatory:

For research use only (RUO)

This Product is for Research Use only (RUO) and not to be used in Diagnostic Procedure.