

## Anti-Denosumab Polyclonal Antibody

Cat# AB3015-100, AB3015-500

### Product description

---

Product: Anti-Denosumab Polyclonal Antibody

Host: Rabbit

Isotype: IgG

Format: Protein-A affinity Purified and depleted against Human IgG.

Immunogen: Denosumab

Product Type: Polyclonal Antibody

Other Names: Xgeva®

### Formulation

---

Anti-Denosumab Polyclonal Antibody stored in 50mM Phosphate Buffer Saline, pH 7.4

### Application

---

The product is purified on protein-A affinity resin and specifically recognizes the denosumab drug. The antibody can be used as positive control for immunogenicity assay. It is recommended that the end user optimizes the product for their particular application with appropriate negative and positive controls.

### Specificity

---

This Anti-Denosumab Polyclonal Antibody recognizes Denosumab specifically.

### Tested Application

---

This product has been tested for ELISA/Western application.

### Storage

---

Store at: -20°C in aliquots and avoid repeated freeze/thaw cycles. If required, 20% sterile Glycerol could be added for long term storage.

### Pack size available

---

100 µg, 500 µg

### Regulatory:

---

For research use only (RUO)

*This Product is for Research Use only (RUO) and not to be used in Diagnostic Procedure.*