

Anti-Insulin Monoclonal Antibody (5G5)

Cat# AB3020-100, AB3020-500

Product description

Product: Anti-Insulin Monoclonal Antibody (5G5)
Host: Mouse
Isotype: IgG
Species: Human
Format: Purified
Immunogen: Insulin
Product Type: Monoclonal Antibody
Other Names: Insulin

Formulation

Anti-Insulin Monoclonal Antibody stored in 1X Phosphate Buffer Saline, 20% Glycerol.

Application

The product is purified on protein-G affinity resin and specifically recognizes the Human Insulin drug. The antibody can be used for bioanalytical assay and western blot. It is recommended that the end user optimizes the product for their particular application with appropriate negative and positive controls.

Specificity

This Anti Insulin antibody recognizes specifically the Human Insulin antibody.

Tested Application

This product has been tested for ELISA and western blot application.

Storage

Store at: -20°C in aliquots and avoid repeated freeze/thaw cycles.

Pack size available

100 µg, 500 µg

Regulatory:

For research use only (RUO)

This Product is for Research Use only (RUO) and not to be used in Diagnostic Procedure.