
Mouse anti Matrix Metalloproteinase 9 (MMP9) Antibody

Cat# AB3075-100

Product description

Product: Mouse anti Matrix Metalloproteinase 9 (MMP9) Antibody
Host: Mouse
Species: Human
Isotype: IgG2b Kappa
Format: Protein A + Protein G Affinity Purified
Immunogen: Recombinant MMP9 (Gly213~Ala399) expressed in E.coli
Product Type: Monoclonal Antibody

Formulation

0.01M PBS, pH7.4, containing 0.05% Proclin-300, 50% glycerol.

Application

The product is affinity purified and specifically recognizes the Human Matrix Metalloproteinase 9 (MMP9). The antibody can be used for variety of application. It is recommended that the end user optimizes the product for their particular application with appropriate controls.

Specificity

The antibody is a mouse monoclonal antibody raised against Matrix Metalloproteinase 9 (MMP9). It has been selected for its ability to recognize Matrix Metalloproteinase 9 (MMP9) in immunohistochemical staining and western blotting.

Tested Application

Western blotting: 0.5-3µg/mL;
Immunohistochemistry: 5-30µg/mL;
Immunocytochemistry: 5-30µg/mL;
Optimal working dilutions must be determined by end user.

Storage

Store at: -20°C in aliquots and avoid repeated freeze/thaw cycles.

Pack size available

100 µg, 500 µg

Regulatory:

For research use only (RUO)

This Product is for Research Use only (RUO) and not to be used in Diagnostic Procedure.

No. A-112, KSSIDC Block-3, Electronic city phase-1, Bangalore 560100, Karnataka, INDIA.

Web: denovobiolabs.com, Email: info@denovobiolabs.com, Contact: +91-80-29575711

Version-1.0
