

---

**Rabbit Anti Solute Carrier Family 30 Member 8 (SLC30A8) Antibody**

Cat# AB3078-100

**Product description**

---

Product: Rabbit Anti Solute Carrier Family 30 Member 8 (SLC30A8) Antibody

Host: Rabbit

Species: Human

Format: Protein A Affinity Purified.

Immunogen: Recombinant SLC30A8 (Thr263~Asp369) expressed in E.coli

Product Type: Polyclonal Antibody

**Formulation**

---

0.01M PBS, pH7.4, containing 0.05% Proclin-300, 50% glycerol.

**Application**

---

The product is affinity purified and specifically recognizes the Human Solute Carrier Family 30 Member 8 (SLC30A8). The antibody can be used for variety of application. It is recommended that the end user optimizes the product for their particular application with appropriate controls.

**Specificity**

---

The antibody is a rabbit polyclonal antibody raised against Solute Carrier Family 30 Member 8 (SLC30A8). It has been selected for its ability to recognize Solute Carrier Family 30 Member 8 (SLC30A8) in immunohistochemical staining and western blotting.

**Tested Application**

---

Western blotting: 0.5-3µg/mL

Immunohistochemistry: 5-30µg/mL

Immunocytochemistry: 5-30µg/mL

Optimal working dilutions must be determined by end user.

**Storage**

---

Store at: -20°C in aliquots and avoid repeated freeze/thaw cycles.

**Pack size available**

---

100 µg, 500 µg

**Regulatory:**

---

For research use only (RUO)

*This Product is for Research Use only (RUO) and not to be used in Diagnostic Procedure.*