
Rabbit Anti Beclin 1 (BECN1) Antibody

Cat# AB3080

Product description

Product: Rabbit Anti Beclin 1 (BECN1) Antibody

Host: Rabbit

Species: Human

Format: Affinity Purified.

Immunogen: Recombinant BECN1 (Met1~Thr273 (Accession # Q14457)) expressed in E.coli

Product Type: Polyclonal Antibody

Formulation

0.01M PBS, pH7.4, containing 0.05% Proclin-300, 50% glycerol.**Application**

The product is affinity purified and specifically recognizes the Human Beclin 1 (BECN1). The antibody can be used for variety of application. It is recommended that the end user optimizes the product for their particular application with appropriate controls.**Specificity**

The antibody is a rabbit polyclonal antibody raised against Beclin 1 (BECN1). It has been selected for its ability to recognize Beclin 1 (BECN1) in immunohistochemical staining and western blotting.**Tested Application**

Western blotting: 0.5-3µg/mL;

Immunohistochemistry: 5-30µg/mL;

Immunocytochemistry: 5-30µg/mL;

Optimal working dilutions must be determined by end user.

Storage

Store at: -20°C in aliquots and avoid repeated freeze/thaw cycles.**Pack size available**

100 µg, 500 µg**Regulatory:**

For research use only (RUO)

This Product is for Research Use only (RUO) and not to be used in Diagnostic Procedure.