

---

**Rabbit Anti Cross Linked C-Telopeptide Of  
Type I Collagen (CTXI) Antibody**  
Cat# AB3084-100

**Product description**

---

Product: Rabbit Anti Cross Linked C-Telopeptide Of Type I Collagen (CTXI) Antibody  
Host: Rabbit  
Species: Human  
Format: Protein A Affinity Purified  
Immunogen: Synthetic Peptide, CTXI conjugated to KLH. Target peptide sequence: EKAHDGGRYRA  
Product Type: Polyclonal Antibody

**Formulation**

---

0.01M PBS, pH7.4, containing 0.05% Proclin-300, 50% glycerol.

**Application**

---

The product is affinity purified and specifically recognizes the Human Cross Linked C-Telopeptide Of Type I Collagen (CTXI). The antibody can be used for variety of application. It is recommended that the end user optimizes the product for their particular application with appropriate controls.

**Specificity**

---

The antibody is a rabbit polyclonal antibody raised against Cross Linked C-Telopeptide Of Type I Collagen (CTXI). It has been selected for its ability to recognize Cross Linked C-Telopeptide Of Type I Collagen (CTXI) in immunohistochemical staining and western blotting.

**Tested Application**

---

Western blotting: 0.5-2µg/mL;  
Immunohistochemistry: 5-20µg/mL;  
Immunocytochemistry: 5-20µg/mL;  
Optimal working dilutions must be determined by end user.

**Storage**

---

Store at: -20°C in aliquots and avoid repeated freeze/thaw cycles.

**Pack size available**

---

100 µg, 500 µg

**Regulatory:**

---

For research use only (RUO)

*This Product is for Research Use only (RUO) and not to be used in Diagnostic Procedure.*