

---

## **Rabbit Anti Cyclin D1 (CCND1) Antibody**

Cat# AB3092-100

### **Product description**

---

Product: Rabbit Anti Cyclin D1 (CCND1) Antibody.  
Host: Rabbit  
Species: Human  
Format: Protein A Affinity Purified  
Immunogen: Recombinant CCND1 (Met1~Ile295) expressed in E.coli.  
Product Type: Polyclonal Antibody

### **Formulation**

---

0.01M PBS, pH7.4, containing 0.05% Proclin-300, 50% glycerol.

### **Application**

---

The product is affinity purified and specifically recognizes the Human Cyclin D1 (CCND1). The antibody can be used for variety of application. It is recommended that the end user optimizes the product for their particular application with appropriate controls.

### **Specificity**

---

The antibody is a rabbit polyclonal antibody raised against Cyclin D1 (CCND1). It has been selected for its ability to recognize CCND1 in immunohistochemical staining and western blotting.

### **Tested Application**

---

Western blotting: 0.5-2µg/mL;  
Immunohistochemistry: 5-20µg/mL;  
Immunocytochemistry: 5-20µg/mL;  
Optimal working dilutions must be determined by end user.

### **Storage**

---

Store at: -20°C in aliquots and avoid repeated freeze/thaw cycles.

### **Pack size available**

---

100 µg, 500 µg

### **Regulatory:**

---

For research use only (RUO)

*This Product is for Research Use only (RUO) and not to be used in Diagnostic Procedure.*