

---

**Mouse Anti Ki-67 Protein (Ki67P) Antibody**

Cat# AB3096-100

**Product description**

---

Product: Mouse Anti Ki-67 Protein (Ki67P) Antibody  
Host: Mouse  
Species: Human  
Isotype: IgG2b Kappa  
Format: Protein A/G Affinity Purified  
Immunogen: RPC047Hu02-Recombinant Ki-67 Protein (Ki67P).  
Product Type: Monoclonal Antibody

**Formulation**

---

0.01M PBS, pH7.4, containing 0.05% Proclin-300, 50% glycerol.

**Application**

---

The product is affinity purified and specifically recognizes the Human Ki-67 Protein (Ki67P). The antibody can be used for variety of application. It is recommended that the end user optimizes the product for their particular application with appropriate controls.

**Specificity**

---

The antibody is a mouse monoclonal antibody raised against Ki-67 Protein (Ki67P). It has been selected for its ability to recognize Ki67P in immunohistochemical staining and western blotting.

**Tested Application**

---

Western blotting: 0.5-5ug/ml  
Immunocytochemistry in formalin fixed cells: 5-30ug/ml  
Immunohistochemistry in formalin fixed frozen section: 5-30ug/ml  
Immunohistochemistry in paraffin section: 5-30ug/ml  
Enzyme-linked Immunosorbent Assay: 0.05-2ug/ml  
Optimal working dilutions must be determined by end user.

**Storage**

---

Store at: -20°C in aliquots and avoid repeated freeze/thaw cycles.

**Pack size available**

---

100 µg, 500 µg

**Regulatory:**

---

For research use only (RUO)

*This Product is for Research Use only (RUO) and not to be used in Diagnostic Procedure.*