
Rabbit Anti Cyclin Dependent Kinase Inhibitor 1A (CDKN1A) Antibody

Cat# AB3102-100

Product description

Product: Rabbit Anti Cyclin Dependent Kinase Inhibitor 1A (CDKN1A) Antibody
Host: Rabbit
Species: Human
Format: Protein A Affinity Purified
Immunogen: Recombinant CDKN1A (Met1~Pro164) expressed in E.coli.
Product Type: Polyclonal Antibody

Formulation

0.01M PBS, pH7.4, containing 0.05% Proclin-300, 50% glycerol.

Application

The product is affinity purified and specifically recognizes the Human Cyclin Dependent Kinase Inhibitor 1A (CDKN1A). The antibody can be used for variety of application. It is recommended that the end user optimizes the product for their particular application with appropriate controls.

Specificity

The antibody is a rabbit polyclonal antibody raised against Cyclin Dependent Kinase Inhibitor 1A (CDKN1A). It has been selected for its ability to recognize CDKN1A in immunohistochemical staining and western blotting.

Tested Application

Western blotting: 0.5-3µg/mL;
Immunohistochemistry: 5-30µg/mL;
Immunocytochemistry: 5-30µg/mL;
Optimal working dilutions must be determined by end user.

Storage

Store at: -20°C in aliquots and avoid repeated freeze/thaw cycles.

Pack size available

100 µg, 500 µg

Regulatory:

For research use only (RUO)

This Product is for Research Use only (RUO) and not to be used in Diagnostic Procedure.