
Rabbit Anti Glypican 1 (GPC1) Antibody

Cat# AB3112-100

Product description

Product: Rabbit Anti Glypican 1 (GPC1) Antibody
Host: Rabbit
Species: Human
Format: Protein A Affinity purified.
Immunogen: Recombinant GPC1 (Leu114~Lys277) expressed in E.coli
Product Type: Polyclonal Antibody

Formulation

PBS, pH7.4, containing 0.02% NaN₃, 50% glycerol.

Application

The product is affinity purified and specifically recognizes the Human Glypican 1 (GPC1). The antibody can be used for variety of application. It is recommended that the end user optimizes the product for their particular application with appropriate controls.

Specificity

The antibody is rabbit polyclonal antibody raised against Glypican 1 (GPC1). It has been selected for its ability to recognize GPC1 in immunohistochemical staining and western blotting.

Tested Application

Western blotting: 0.5-2µg/mL;
Immunohistochemistry: 5-20µg/mL;
Immunocytochemistry: 5-20µg/mL;
Optimal working dilutions must be determined by end user.

Storage

Store at: -20°C in aliquots and avoid repeated freeze/thaw cycles.

Pack size available

100 µg, 500 µg

Regulatory:

For research use only (RUO)

This Product is for Research Use only (RUO) and not to be used in Diagnostic Procedure.