
Mouse Anti Hyaluronic Acid (HA) Antibody

Cat# AB3120-100

Product description

Product: Mouse Anti Hyaluronic Acid (HA) Antibody

Host: Mouse

Species: Pan species (General)

Isotype: IgG1 Kappa

Format: Protein A + Protein G Affinity purified.

Immunogen: Small Molecule, HA conjugated to OVA.

Product Type: Monoclonal Antibody

Formulation

PBS, pH7.4, containing 0.02% NaN₃, 50% glycerol.**Application**

The product is affinity purified and specifically recognizes the Hyaluronic Acid (HA). The antibody can be used for variety of application. It is recommended that the end user optimizes the product for their particular application with appropriate controls.

Specificity

The antibody is mouse monoclonal antibody raised against Hyaluronic Acid (HA). It has been selected for its ability to recognize HA in immunohistochemical staining and western blotting.

Tested Application

Immunohistochemistry: 5-20µg/mL;

Immunocytochemistry: 5-20µg/mL;

Optimal working dilutions must be determined by end user.

Storage

Store at: -20°C in aliquots and avoid repeated freeze/thaw cycles.

Pack size available

100 µg, 500 µg

Regulatory:

For research use only (RUO)

This Product is for Research Use only (RUO) and not to be used in Diagnostic Procedure.