
Mouse Anti Hexa Histidine (H6) Antibody

Cat# AB3123-100

Product description

Product: Mouse Anti Hexa Histidine (H6) Antibody
Host: Mouse
Species: Pan species (General)
Isotype: IgG1 Kappa
Format: Protein A + Protein G Affinity purified.
Immunogen: Small Molecule, H6 conjugated to BSA.
Product Type: Monoclonal Antibody

Formulation

0.01M PBS, pH7.4, containing 0.05% Proclin-300, 50% glycerol.

Application

The product is affinity purified and specifically recognizes the Hexa Histidine (H6). The antibody can be used for variety of application. It is recommended that the end user optimizes the product for their particular application with appropriate controls.

Specificity

The antibody is mouse monoclonal antibody raised against Hexa Histidine (H6). It has been selected for its ability to recognize H6 in immunohistochemical staining and western blotting.

Tested Application

Immunohistochemistry: 5-20µg/mL;
Immunocytochemistry: 5-20µg/mL;
Optimal working dilutions must be determined by end user.

Storage

Store at: -20°C in aliquots and avoid repeated freeze/thaw cycles.

Pack size available

100 µg, 500 µg

Regulatory:

For research use only (RUO)

This Product is for Research Use only (RUO) and not to be used in Diagnostic Procedure.