
Mouse Anti Lipopolysaccharide (LPS) Antibody

Cat# AB3124-100

Product description

Product: Mouse Anti Lipopolysaccharide (LPS) Antibody

Host: Mouse

Species: Pan species (General)

Isotype: IgG2b Kappa

Format: Protein A + Protein G Affinity purified.

Immunogen: Small Molecule, LPS conjugated to OVA.

Product Type: Monoclonal Antibody

Formulation

0.01M PBS, pH7.4, containing 0.05% Proclin-300, 50% glycerol.

Application

The product is affinity purified and specifically recognizes the Lipopolysaccharide (LPS). The antibody can be used for variety of application. It is recommended that the end user optimizes the product for their particular application with appropriate controls.

Specificity

The antibody is mouse monoclonal antibody raised against Lipopolysaccharide (LPS). It has been selected for its ability to recognize LPS in immunohistochemical staining and western blotting.

Tested Application

Immunohistochemistry: 5-20µg/mL;

Immunocytochemistry: 5-20µg/mL;

Optimal working dilutions must be determined by end user.

Storage

Store at: -20°C in aliquots and avoid repeated freeze/thaw cycles.

Pack size available

100 µg, 500 µg

Regulatory:

For research use only (RUO)

This Product is for Research Use only (RUO) and not to be used in Diagnostic Procedure.