

---

**Rabbit Anti Putrescine (PT) Antibody**

Cat# AB3126-100

**Product description**

---

Product: Rabbit Anti Putrescine (PT) Antibody  
Host: Rabbit  
Species: Pan species (General)  
Format: Protein A Affinity purified.  
Immunogen: Small Molecule, Putrescine conjugated to OVA.  
Product Type: Polyclonal Antibody

**Formulation**

---

0.01M PBS, pH7.4, containing 0.05% Proclin-300, 50% glycerol.

**Application**

---

The product is affinity purified and specifically recognizes the Putrescine (PT). The antibody can be used for variety of application. It is recommended that the end user optimizes the product for their particular application with appropriate controls.

**Specificity**

---

The antibody is rabbit polyclonal antibody raised against Putrescine (PT). It has been selected for its ability to recognize PT in immunohistochemical staining and western blotting.

**Tested Application**

---

Enzyme-Linked Immune Sorbent Assay: 100ng/mL;  
Immunohistochemistry: 5-20µg/mL;  
Immunocytochemistry: 5-20µg/mL;  
Optimal working dilutions must be determined by end user.

**Storage**

---

Store at: -20°C in aliquots and avoid repeated freeze/thaw cycles.

**Pack size available**

---

100 µg, 500 µg

**Regulatory:**

---

For research use only (RUO)

*This Product is for Research Use only (RUO) and not to be used in Diagnostic Procedure.*