

---

**Rabbit Anti Oxytocin (OT) Antibody**

Cat# AB3129-100

**Product description**

---

Product: Rabbit Anti Oxytocin (OT) Antibody  
Host: Rabbit  
Species: Pan species (General)  
Format: Protein A Affinity purified.  
Immunogen: Small Molecule, Oxytocin conjugated to OVA.  
Product Type: Polyclonal Antibody

**Formulation**

---

0.01M PBS, pH7.4, containing 0.05% Proclin-300, 50% glycerol.

**Application**

---

The product is affinity purified and specifically recognizes the Oxytocin (OT). The antibody can be used for variety of application. It is recommended that the end user optimizes the product for their particular application with appropriate controls.

**Specificity**

---

The antibody is rabbit polyclonal antibody raised against Oxytocin (OT). It has been selected for its ability to recognize OT in immunohistochemical staining and western blotting.

**Tested Application**

---

Immunohistochemistry: 5-20µg/mL;  
Immunocytochemistry: 5-20µg/mL;  
Optimal working dilutions must be determined by end user.

**Storage**

---

Store at: -20°C in aliquots and avoid repeated freeze/thaw cycles.

**Pack size available**

---

100 µg, 500 µg

**Regulatory:**

---

For research use only (RUO)

*This Product is for Research Use only (RUO) and not to be used in Diagnostic Procedure.*