

---

**Rabbit Anti Bilirubin (Bb) Antibody**

Cat# AB3133-100

**Product description**

---

Product: Rabbit Anti Bilirubin (Bb) Antibody  
Host: Rabbit  
Species: Pan species (General)  
Format: Protein A Affinity purified.  
Immunogen: Small Molecule, Bilirubin conjugated to OVA.  
Product Type: Polyclonal Antibody

**Formulation**

---

0.01M PBS, pH7.4, containing 0.05% Proclin-300, 50% glycerol.

**Application**

---

The product is affinity purified and specifically recognizes the Histamine (HA). The antibody can be used for variety of application. It is recommended that the end user optimizes the product for their particular application with appropriate controls.

**Specificity**

---

The antibody is rabbit polyclonal antibody raised against Histamine (HA). It has been selected for its ability to recognize HA in immunohistochemical staining and western blotting.

**Tested Application**

---

Enzyme-Linked Immune Sorbent Assay: 4.53µg/mL;  
Immunohistochemistry: 5-20µg/mL;  
Immunocytochemistry: 5-20µg/mL;  
Optimal working dilutions must be determined by end user.

**Storage**

---

Store at: -20°C in aliquots and avoid repeated freeze/thaw cycles.

**Pack size available**

---

100 µg, 500 µg

**Regulatory:**

---

For research use only (RUO)

---

*This Product is for Research Use only (RUO) and not to be used in Diagnostic Procedure.*

No. A-112, KSSIDC Block-3, Electronic city phase-1, Bangalore 560100, Karnataka, INDIA.

Web: [denovobiolabs.com](http://denovobiolabs.com), Email: [info@denovobiolabs.com](mailto:info@denovobiolabs.com), Contact: +91-80-29575711

Version-1.0

---