
Rabbit Anti Vitamin B2 (VB2) Antibody

Cat# AB3134-100

Product description

Product: Rabbit Anti Vitamin B2 (VB2) Antibody
Host: Rabbit
Species: Pan species (General)
Isotype: IgG
Format: Affinity purified.
Immunogen: VB2- OVA.
Product Type: Polyclonal Antibody

Formulation

PBS, pH7.4, containing 0.02% NaN₃, 50% glycerol.

Application

The product is affinity purified and specifically recognizes the Vitamin B2 (VB2). The antibody can be used for variety of application. It is recommended that the end user optimizes the product for their particular application with appropriate controls.

Specificity

The antibody is rabbit polyclonal antibody raised against Vitamin B2 (VB2). It has been selected for its ability to recognize estrogen receptor in immunohistochemical staining and western blotting.

Tested Application

Western blotting: 1:100-400
Immunocytochemistry in formalin fixed cells: 1:100-500
Immunohistochemistry in formalin fixed frozen section: 1:100-500
Immunohistochemistry in paraffin section: 1:50-200
Enzyme-linked Immunosorbent Assay: 1:100-200
Optimal working dilutions must be determined by end user.

Storage

Store at: -20°C in aliquots and avoid repeated freeze/thaw cycles.

Pack size available

100 µg, 500 µg

Regulatory:

For research use only (RUO)

This Product is for Research Use only (RUO) and not to be used in Diagnostic Procedure.