
Rabbit Anti Progesterone Receptor Antibody

Cat# AB3137-100

Product description

Product: Rabbit Anti Progesterone Receptor Antibody
Host: Rabbit
Species: Pan species(General)
Isotype: IgG
Format: Affinity Purified
Immunogen: Recombinant Human Progesterone Receptor
Product Type: Polyclonal Antibody

Formulation

PBS, pH7.4, containing 0.02% NaN₃, 50% glycerol.

Application

The product is affinity purified and specifically recognizes the Human Progesterone Receptor. The antibody can be used for variety of application. It is recommended that the end user optimizes the product for their particular application with appropriate controls.

Specificity

The antibody is a rabbit polyclonal antibody raised against PGR. It has been selected for its ability to recognize PGR in immunohistochemical staining and western blotting.

Tested Application

Western blotting: 0.5-2µg/mL;
Immunohistochemistry: 5-20µg/mL;
Immunocytochemistry: 5-20µg/mL;
Optimal working dilutions must be determined by end user.

Storage

Store at: -20°C in aliquots and avoid repeated freeze/thaw cycles.

Pack size available

100 µg, 500 µg

Regulatory:

For research use only (RUO)

This Product is for Research Use only (RUO) and not to be used in Diagnostic Procedure.