

Goat Anti Human IgG Fc-specific-HRP Antibody

Cat# AC001

Product description

Product: Goat Anti Human IgG Fc-specific-HRP Antibody

Host: Goat

Isotype: IgG

Format: Affinity Purified

Product Type: Polyclonal Antibody Conjugate

Other Names: NA

Formulation

Goat Anti Human IgG Fc-specific-HRP Antibody stored in 10 mM Phosphate Buffer Saline pH 7.4, 10 mg/ml BSA, 50% Glycerol.

Application

The product is purified on affinity resin and specifically recognizes the Fc fragment of human IgG and tested for ELISA and western blot application. It is recommended that the end user optimizes the product for their particular application with appropriate negative and positive controls.

Specificity

This Goat Anti Human IgG Fc-specific-HRP Antibody recognizes specifically the Fc fragment of human IgG but not the Fab fragment of human IgG. The antibody may cross react with immunoglobulins from other species.

Tested Application

Suggested dilution for ELISA 1: 2000 - 1:10000 and western blot 1: 2000 - 1:5000.

Storage

Store at: -20°C in aliquots and avoid repeated freeze/thaw cycles.

Pack size available

1 ml

Regulatory:

For research use only (RUO)

This Product is for Research Use only (RUO) and not to be used in Diagnostic Procedure.