
Goat Anti Rabbit-HRP Antibody

Cat# AC003

Product description

Product: Goat Anti Rabbit-HRP Antibody
Host: Goat
Isotype: IgG
Format: Affinity Purified
Product Type: Polyclonal Antibody Conjugate
Other Names: NA

Formulation

Goat Anti Rabbit-HRP Antibody stored in 10 mM Phosphate Buffer Saline pH 7.4, 10 mg/ml BSA, Conjugate Stabilizer.

Application

The product is purified on affinity resin and specifically recognizes the Rabbit IgG and tested for ELISA and western blot application. It is recommended that the end user optimizes the product for their particular application with appropriate negative and positive controls.

Specificity

This Goat Anti Rabbit-HRP Antibody recognizes specifically the Rabbit IgG. The antibody may cross react with immunoglobulins from other species.

Tested Application

Suggested dilution; for ELISA 1: 1000 - 1:5000 and Western blot 1: 500 - 1:2000.

Storage

Store at: -20°C in aliquots and avoid repeated freeze/thaw cycles.

Pack size available

1 ml

Regulatory:

For research use only (RUO)

This Product is for Research Use only (RUO) and not to be used in Diagnostic Procedure.