

---

**Goat Anti Guinea Pig-HRP Antibody**

Cat# AC004

**Product description**

---

Product: Goat Anti Guinea Pig-HRP Antibody

Host: Goat

Isotype: IgG

Format: Affinity Purified

Product Type: Polyclonal Antibody Conjugate

Other Names: NA

**Formulation**

---

Goat Anti Guinea Pig-HRP Antibody stored in 10 mM Phosphate Buffer Saline pH 7.4, 10 mg/ml BSA, 50% Glycerol.

**Application**

---

The product is purified on affinity resin and specifically recognizes the Guinea Pig IgG and tested for ELISA and western blot application. It is recommended that the end user optimizes the product for their particular application with appropriate negative and positive controls.

**Specificity**

---

This Goat Anti Guinea Pig -HRP Antibody recognizes specifically the Guinea Pig IgG. The antibody may cross react with immunoglobulins from other species.

**Tested Application**

---

Suggested dilution; for ELISA 1: 1000 - 1:5000 and Western blot 1: 500 - 1:2000.

**Storage**

---

Store at: -20°C in aliquots and avoid repeated freeze/thaw cycles.

**Pack size available**

---

1 ml

**Regulatory:**

---

For research use only (RUO)

*This Product is for Research Use only (RUO) and not to be used in Diagnostic Procedure.*