
Rabbit Anti Goat-HRP Antibody
Cat# AC005**Product description**

Product: Rabbit Anti Goat-HRP Antibody
Host: Rabbit
Isotype: IgG
Format: Affinity Purified
Product Type: Polyclonal Antibody Conjugate
Other Names: NA

Formulation

Rabbit Anti Goat-HRP Antibody stored in 10 mM Phosphate Buffer Saline pH 7.4, 10 mg/ml BSA, 50% Glycerol.

Application

The product is purified on affinity resin and specifically recognizes the Goat IgG and tested for ELISA and western blot application. It is recommended that the end user optimizes the product for their particular application with appropriate negative and positive controls.

Specificity

This Rabbit Anti Goat-HRP Antibody recognizes specifically the Goat IgG. The antibody may cross react with immunoglobulins from other species.

Tested Application

Suggested dilution; for ELISA 1: 1000 - 1:5000 and Western blot 1: 500 - 1:2000.

Storage

Store at: -20°C in aliquots and avoid repeated freeze/thaw cycles.

Pack size available

1 ml

Regulatory:

For research use only (RUO)

This Product is for Research Use only (RUO) and not to be used in Diagnostic Procedure.