

---

**Streptavidin-Peroxidase Conjugate**

Cat# AC006-1ML

**Product description**

---

Product: Streptavidin- peroxidase conjugate  
Pack Size: 1 ml  
Form: Liquid

**Formulation**

---

0.01M Sodium Phosphate, 0.25M NaCl, pH 7.6, 10 mg/ml BSA, Conjugate Stabilizer.

**Application**

---

The product is tested for ELISA, western blot with ECL, immunocytochemistry application. It is recommended that the end user optimizes the product for their particular application with appropriate negative and positive controls.

**Tested Application**

---

Suggested dilution; for ELISA 1: 500 - 1:1000 and Western blot 1: 250 - 1:500.

**Storage**

---

Store at: -20°C in aliquots and avoid repeated freeze/thaw cycles.

**Pack size available**

---

1 ml

**Regulatory:**

---

For research use only (RUO)

*This Product is for Research Use only (RUO) and not to be used in Diagnostic Procedure.*