

## PRODUCT INFORMATION

**Product** : PBS Solution [1X] Ready To Use, pH 7.4  
**Catalog No.** : AR1015-100ML  
**Form** : Liquid  
**Pack Size** : 100 ml , 500 ml  
**Storage** : 2-8 °C  
**Shipping Temperature:** Ambient temperature

## PRODUCT DESCRIPTION

Phosphate Buffered Saline (PBS) Solution is widely used in biological research. The reagent is Clear Solution. It contains monobasic potassium phosphate, sodium chloride, potassium chloride, and dibasic sodium phosphate. This solution does not contain any preservative components.

## DIRECTION FOR USE

- Buffer is a clear solution.
- Ready to use.

## PRECAUTIONS

- Do not dilute the solution
- Avoid inter lot/batch mixing

## REGULATORY

For research use only (RUO)

*The end user must determine the suitability of the product for its particular use. For any technical query, mail us at [info@denovobiolabs.com](mailto:info@denovobiolabs.com);  
**Research Purpose Only***